

ViewPoint

Connect. Report. Relax.

Data Sheet Scanners with Transfer of Measurement Data



In this data sheet all ultrasound machines are listed that are capable of transferring measurement data to ViewPoint software.

It is the responsibility of the customer to make sure that the measurement data transfer option is enabled on the ultrasound machine.

Serial

GE Medical Systems

LogiqBook
Logiq 3
Logiq 5 (BT04 and lower)
Logiq 7
Logiq 9
Logiq 200 ProW
Logiq 400
Logiq 500
Logiq 500 Pro
Logiq 700 ES
Logiq 700 MR

Aloka

SSD-1000
SSD-1400
SSD-4000
SSD-5000
SSD-5500

Hitachi

EUB 525

Philips Medical Systems

Envisor
HD11
HDI 3000
HDI 4000
HDI 5000
Ultramark 9

Siemens

Antares
G 50
G 60
Sonoline Elegra
Sonoline Omnia
Sonoline Prima
Sonoline Sienna
Sonoline Versa
128XP
Aspen
Sequoia 512

Sonoace/Medison

Accuvix XQ
SA8000Live

Toshiba

CAPASEE SSA-220A
CoreVision SSH-350A
EcoCee SSH-340A
JustVision 200 SSA-320A
JustVision 400 SSA-325A
Nemio
PowerVision 6000 SSA-370A
PowerVision 7000 SSA-380A
PowerVision 8000
SSH-140 A
SSH-160A
SSH-270A
SSA-250A
SSA-270A
TOSBEE SSA-340A

DICOM SR

GE Medical Systems

Logiq P5 (BT06)
Logiq 9 (BT04 and higher)
Voluson 730 (BT05 and higher)
Voluson i/Volsuon e (BT06 and higher)
Voluson E8 (BT06 and higher)
Vivid 7 (BT04 and higher)

Philips Medical Systems

HD11XE
iU22

File Transfer

GE Medical Systems

Vivid 3
Vivid 4
Vivid 7(BT03 and lower)
Vivid i
Voluson 530 D
Voluson 730

Hitachi

EUB 6000
EUB 6500

Toshiba

Xario
Aplio

